SEVERIN

Gebrauchsanweisung	Fußmassagegerät (D)
Instructions for use	Foot massager (GB)
Mode d'emploi	Thalasso-pieds (F)
Gebruiksaanwijzing	Voetmassage apparaat (NL)
Instrucciones de uso	Masajeador de pies electrónico (E)
Manuale d'uso	Idromassaggiatore per piedi (I)
Brugsanvisning	Fodmassageapparat (DK)
Bruksanvisning	Fotbad med massage (S)
Käyttöohje	Jalkahierontalaite (FIN)
Instrukcja obsługi	Aparat do masażu stóp (PL)
Οδηγίες χρήσεως	Συσκευή για μασάζ ποδιών (GR)
Руководство по эксплуатации	Ванна для массажа ног (RUS)

Lieber Kunde!

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Foot massager

Dear Customer,

Before using the appliance, the user should read the following instructions carefully.

Connection to the mains supply

Make sure that the supply voltage corresponds with the voltage marked on the rating label of the appliance. This product complies with all binding CE labelling directives.

Familiarisation

1. Storage compartment for pedicure unit
2. Rotary control
3. Cover
4. Nail-polish dryer
5. Foot massager
6. Roller massage insert for the soles
7. Infrared pad
8. Water-level indicator
9. Right foot-pad with massage studs
10. Power cord with plug
11. Emptying recess
12. Slide switch
13. Battery compartment cover
14. Pedicure unit
15. Aroma chamber
16. Sponge attachment
17. Abrading attachment
18. Filing attachment
19. Brush attachment
20. Left foot-pad with massage studs
21. Rating label (on the side)

Important safety instructions

- Before the appliance is used, the main body as well as any attachments fitted should be checked thoroughly for any defects. If the appliance, for instance, has been dropped onto a hard surface, it must no longer be used: even invisible damage may have adverse effects on its operational safety.
- If the appliance is being used in a

bathroom, it must be positioned outside the 0.6 m area (see diagram) in such a way that the rotary control cannot be reached by anyone using the shower or bath.

- Additional protection is provided by fitting an earth-leakage circuit-breaker, with a rated tripping current not exceeding 30mA, into your bathroom installation. Please consult a qualified electrician.
- Do not attempt to charge or by any other means reactivate the batteries of the pedicure unit; never take batteries apart, or burn or short-circuit them.
- Always remove the batteries from the pedicure unit once they are exhausted, and also during extended periods of non-use. This will prevent damage from leakage.
- All batteries in the pedicure unit should be replaced at the same time. Do not use batteries of different types or brands, or those of different capacities.
- Do not leave the foot massager unattended while it is switched on, nor let children use it without supervision.
- Do not stand on the foot massager; use it only while sitting.
- The foot massager becomes hot during use. Those who are over-sensitive to heat should take extreme care when using the appliance.
- Do not operate the foot massager near a sink, shower or bath.
- This appliance is intended for domestic use only and is not designed for medical or commercial applications.
- Always switch the foot massager off and remove the plug from the wall socket after use and during cleaning, and also in case of any malfunction. When removing the plug from the wall socket, ensure that your hands are dry.
- The power cord should be regularly examined for any signs of damage. In the event of such damage being found, the appliance must no longer be used.

- Always switch the pedicure unit off before changing any attachments.
- This appliance should only be used with the attachments provided by the manufacturer.
- Do not use the pedicure unit if you suffer from eczema or open wounds on the feet.
- This appliance is not intended for use by any person (including children) with reduced physical, sensory or mental capabilities, or lacking experience and knowledge, unless they have been given supervision or instruction concerning the use of the appliance by a person responsible for their safety.
- Children should be supervised to ensure that they do not play with the appliance.
- Keep children well away from packaging materials, as they are a potential source of danger e.g. from suffocation.
- No responsibility will be accepted if damage results from improper use, or if these instructions are not complied with.
- In order to comply with safety regulations and to avoid hazards, repairs of electrical appliances must be carried out by qualified personnel, including the replacement of the power cord. If repairs are needed, please send the appliance to one of our customer service departments. The address can be found in the appendix to this manual.

Foot massager

General information

A foot massage helps to reduce tension in muscle and tissue areas, enhancing the blood circulation and simultaneously stimulating the metabolic process, thus cleansing the body tissue of waste residues.

- If required, the roller massage insert can be fitted to the foot-pads. Roller massage is specifically recommended for the sole area of the feet.
An infrared pad for heat treatment is located between the foot pads.
- In order to enhance the massage effect, a

foot massage should be applied using warm water. A non-foaming bath gel may also be added.

- The integral heating element located underneath the studded foot-pads prevents the warm water from cooling down too fast, although it is **not** designed to heat up the water.
- The user is recommended to refrain from massage treatment when suffering from phlebitis in the legs, varicose veins (if they tend to bleed easily), skin rash, inflamed or swollen feet or any pain of unknown origin. If in doubt, consult your doctor before using the appliance.

Control settings

- 0** Unit is switched off
- 1** Bubble function + massage + infrared
- 2** Warm air stream for nail dryer
- 3** Massage + infrared

Operation

- Place the foot massager on the floor in front of a chair.
- If required, the massage inserts for the foot-reflex zones can be fitted.
- Fill with warm water until the indicated water level is reached.
- Insert the plug into the wall socket. Make sure that the rotary control is in the “**0**” position.
- Seat yourself comfortably on the chair and position your feet on the foot pads. **Do not** stand on the unit.
- Select the desired function with the rotary control. The recommended application time is 10 to 20 minutes.
- After use, set the control switch to “**0**” again before removing the plug from the wall socket.
- Empty out the water afterwards through the recess.

Nail-polish dryer with aroma function

- Turn the foot massager on and set it to position **2**.
- After applying the nail polish to your toe

- nails, place your feet on the middle pad.
- If the aroma function is required, pull out the nail-dryer insert from the main unit. Add a few drops of scented oil to the sponge inside the aroma chamber, then slide the nail-dryer insert back into the main unit until it locks into place.
- Switch off the foot massager after use.

Pedicure unit

Batteries

This pedicure unit requires 2 type AA / 1.5 V batteries; the use of alkaline batteries is preferable, as they last longer than standard carbon-zinc batteries. The batteries are inserted as follows:

- Slide the battery compartment cover in the direction indicated by the arrow and take it off.
- Insert the batteries according to the schematic diagram inside the compartment, ensuring correct polarity (+/-).
- Replace the compartment cover and slide it back until it is heard locking into place.
- A significant decrease in operating speed is a clear indication that the batteries are exhausted, and that they should be replaced.
- Always remove the batteries from the pedicure unit during extended periods of non-use. This will prevent damage to the appliance from leakage.
- Never dispose of old batteries in your domestic garbage. Old batteries must be disposed of through local ‘special garbage collection points’.

Using the pedicure unit

With the appliance switched off, fit the attachment required onto the front of the pedicure unit and push it in as far as it will go.

- **Sponge attachment:** for polishing toe nails.
- **Abrading attachment:** for shortening and

shaping toe nails.

- **Cleaning brush:** for cleaning the nails after filing, and for removing skin particles remaining after the treatment. It is also useful for preparing the nail surfaces for polishing.
- **Filing attachment:** for filing and other treatment of the nail surfaces.

The slide-switch on the pedicure unit allows for two different speed settings.

The feet should be treated after a hot bath or shower, as calluses and nails will be softened in this way.

Always switch the pedicure unit off after use. A decrease in operating speed indicates low battery capacity. When this occurs, switch off the pedicure unit and replace the batteries.

Cleaning and care

- Before cleaning the appliance, ensure it is disconnected from the power supply.
- To avoid the risk of electric shock, do not clean the housing of the foot massager with water and do not immerse it in water.
- Do not use abrasives or harsh cleaning solutions.

Foot massager

- The inside and the exterior of the foot massager may be cleaned with a slightly damp lint-free cloth. Use household tissue paper to wipe dry afterwards. Special, commercially available cleaning agents for synthetic materials may also be used for cleaning.

Pedicure unit

- The exterior of the pedicure unit may be wiped with a dry cloth.
- The inserts may be cleaned as necessary with an alcohol solution.

Disposal

Do not dispose of old or defective appliances in domestic garbage; this should only be done through public collection points.

Guarantee

This product is guaranteed against defects in materials and workmanship for a period of two years from the date of purchase. Under this guarantee the manufacturer undertakes to repair or replace any parts found to be defective, providing the product is returned to one of our authorised service centres. This guarantee is only valid if the appliance has been used in accordance with the instructions, and provided that it has not been modified, repaired or interfered with by any unauthorised person, or damaged through misuse.

This guarantee naturally does not cover wear and tear, nor breakables such as glass and ceramic items, bulbs etc. If the product fails to operate and needs to be returned, pack it carefully, enclosing your name and address and the reason for return. If within the guarantee period, please also provide the guarantee card and proof of purchase.

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Oy Harry Marcell Ab